



LOGIQ e R8

NextGen High Definition Portable Veterinary Ultrasound



High Resolution HD Imaging

Lightweight & Portable

Automatic Image Optimization

On-Board Veterinary Presets



LOGIQ™ e R8 NextGen is the latest edition of GE's new high performance R-series portable ultrasounds. The NextGen operates on a 128-channel processing system, with a solid state hard drive, allowing for unprecedented imaging resolution and performance in a compact console. The S-agile acoustic ultrasound software platform is adopted from GE's flagship S-series ultrasound series, and combines SRI speckle reduction with CrossXBeam compound imaging for exceptional resolution without sacrificing frame rates. The agile acoustic architecture design uses multiple automated imaging algorithms to reduce the requirement of user optimization of the ultrasound for outstanding imaging.

Features:

Automatic Image Optimization

One-click button to optimize image quality, contrast, colour flow and doppler information.

CrossXBeam and SRI Technology

Provides exceptional image resolution and speckle reduction without negative effects on frame rate.

On-Board Veterinary Presets

For simplified "point and shoot" ultrasound.

Harmonic Imaging, Colour Doppler, and Pulsed Wave Doppler

20 Seconds Boot Up Time

New High-Res PDI Blood Flow Detection

In low velocity settings.

Lightweight Portable Console

Under 10 lbs with 1 hour battery life.

The Only Ultrasound Platform With Advanced Raw Data Capabilities

Capture your images now, and optimize, annotate and measure later!

Available NextGen Software Package

Including HD Directional Power Doppler, Panoramic Imaging (LOGIQview), Needle Detection Software, and the Follow Up Tool.

eR8

Most Common Veterinary Probes

6S-RS Cardiac Phased Array	2.5–7 MHz
L10-22-RS Ultra-High Frequency Small Footprint Linear	10–22 MHz
LK760 Linear Transrectal	3.5–10 MHz
12L-RS Linear	4.2–13 MHz
8C-RS Microconvex	4–10 MHz

Continuing Education

We believe continuing education is the key to success. Our onsite installation and training is performed by Veterinary sonographers, and is tailored to your clinics skill level. Looking to improve your ultrasound skills further? We offer a vast array of online and in-person ultrasound courses including: basic to advanced abdominal ultrasound, organ-focused abdominal ultrasound, basic to advanced echocardiography, F.A.S.T ultrasound and musculoskeletal ultrasound.



We are the exclusive Canadian Veterinary Distributor for GE

www.heskavet.ca +1.866.382.6937

©2022 Heska Corporation. All Rights Reserved. LOGIQ is a trademark of The General Electric Company. HESKA is a trademark of Heska Corporation in the U.S. and other countries. CA20LTO906

